

Glass Stand Only for Plantar Testing

When large group studies, are required IITC offers either version of our glass stands to be purchased individually. Allowing users to quickly and efficiently test large study groups.

The two versions available -6 rats/12 mice 2 rats/4 mice

Both sizes are supplied standard -Black anodized base, 5 black anodized legs. Both sizes are supplied with the standard 8" legs. Tempered glass is placed on top. The clear plexi-glas animal enclosures are placed on top of the glass and allow unobstructed views of animal subjects during testing.



Due to the fact that glass acts as a heat sink; IITC developed the only commercially available heated glass. Warming the glass eliminates the heat sink problem, thus animals are placed onto the slightly warmed glass. Builtin temperature controller in the heated base version keeps glass at a constant temperature between 30-40 deg. C. which can be adjusted by user as needed.

Model 400 Heated Base

Acrylic animal enclosures are supplied, with individual lids for each division, black panels divide the sections of animals.

IITC customizes our animal enclosures when needed to any specified dimensions for a slight additional charge. Please contact IITC with your request and we will be happy to supply you pricing on your custom enclosure needs.

Two year warranty **CE Certified**

Power: 120/240 VAC Switchable

Part #	Description
400	Glass Stand for Plantar testing on 12 mice or 6 rats
400G	Heated Glass Stand for Plantar testing on 12 mice or 6 rats
402	Glass Stand for Plantar testing on 4 mice or 2 rats
402G	Glass Stand for Plantar testing on 12 mice or 6 rats
???	Animal Enclosure